

Configurations possibles sur l'AccuriC6

* Disponible avec le Selectable Lasers Module

Detector	Filter	Fluorochrome
FL1	530/30	FITC, GFP, CFSE
FL2	585/40	PE, PI PE-Texas Red
FL3	670 LP	PerCP-Cy5.5 PE-Cy5, PE-Cy7
FL4	675/25	APC Alexa-647 PE-Cy5

3 blue 1 red: Standard Filters

Detector	Filter	Fluorochrome
FL1	530/30	FITC, GFP, CFSE
FL2	585/40	PE
FL3	610/20	PI PE-Texas Red
FL4	675/25	APC Alexa-647 PE-Cy5

3 blue 1 red: 610/20 in FL3

Detector	Filter	Fluorochrome
FL1	530/30	FITC, GFP, CFSE
FL2	585/40	PE, PI PE-Texas Red
FL3	780/60	APC-Cy7 (and substitutes)
FL4	675/25	APC (and substitutes)

* 2 blue 2 red configuration

Detector	Filter	Fluorochrome
FL1	530/30	FITC, GFP, CFSE
FL2	585/40	PE
FL3	675/25	PerCP-Cy5.5 PE-Cy5
FL4	630/30 610/20	PI PE-Texas Red

* 4 blue: FL3: 675 BP FL4: 630 or 610 BP

Detector	Filter	Fluorochrome
FL1	530/30	FITC, GFP, CFSE
FL2	585/40	PE
FL3	780/60	PE-Cy7 (and substitutes)
FL4	630/30 610/20	PI PE-Texas Red

* 4 blue: FL3: 780 BP FL4: 630 or 610 BP

Detector	Filter	Fluorochrome
FL1	530/30	FITC, GFP, CFSE
FL2	585/40	PE, PI PE-Texas Red
FL3	780/60	PE-Cy7 (and substitutes)
FL4	675/25	PerCP-Cy5.5 PE-Cy5

* 4 blue: FL3: 780 BP FL4: 675 BP